



Shavers, Burs & bipolar RF-Electrodes *designed for UBE/BESS*



StenoBur & StenoVap

REGER Shaver System, Handpiece and Foot Switch



Diamond Bur designed for controlled abrasion and bone smoothing!

Single Use and Reusable

StenoVap							
Description	90° Electrode	70° Electrode	Hook slim Electrode	Hook Electrode	45° Bent Electrode	45° Electrode	45° Bendable Electrode
Electrodes							
Length [mm]	150	150	150	150	150	150	150
Electrodes compatible with common RF-Units!							


Ecologically sustainable



RF-Unit, Handpiece and Foot Switch

StenoBur						
Description	Diamond Round Bur SideGuard	Round Bur SideGuard	Barrel Bur SideGuard	Pit Bull	Wave Cutter	Synovia Cutter
Shaver						
Tube Diameter [mm]	5,5	5,5	5,5	4,0	4,0	4,0
Length [mm]	130	130	130	130	130	130
Shavers & Burs compatible with common manufacturers' systems!						





StenoBur & StenoVap

REGER Medizintechnik GmbH
Gewerbestraße 10
D-78667 Villingendorf
Tel +49 741 270 698-0
www.reger-med.de
info@reger-med.de



REGER
+medizintechnik gmbh

Index B 2023